


## Color measurement of Cosmetics ( Foundation and Eyeshadow)

Color measurement is performed when expressing color quantitatively. This makes it possible to evaluate the brightness of a sample and slight color difference.

We picked up foundation and eye shadow (cosmetics), and evaluated them using U-3900 spectrophotometer,  $\phi$  150 integrating sphere and option package program. By using the color value on the coordinates of the  $L^*a^*b^*$  colorimetric system (by using D65 light source and 2 degree visual field), it is possible to easily check slight differences in brightness and color appearance of foundation.

Samples:



Powder foundation



Eye shadow

Measurement conditions:


Scan Speed: 300 nm/min

Slit: 5 nm


Sampling interval: 1 nm

Standard Reflectance Material: Aluminum Oxide White plate

Measurement method: n-d( JIS Z8722 (2000) ), [0/d] (CIE No.15.2)



Foundation Diffuse Reflectance spectrum



Eye shadow Diffuse Reflectance spectrum

## Color Measurement of Foundations

Color calculation results (D<sub>65</sub> light source, 2 degree visual field)

No.	Sample name	L*	a*	b*
1	Foundation 1	76.73	10.71	22.20
2	Foundation 2	71.89	12.08	21.92
3	Foundation 3	69.11	11.66	22.09

Among the three foundations, Foundation 1 had the highest value of L\* (brightness). From this we can see that Foundation 1 has the strongest brightness.

Foundation 2 has the highest a\* value.  
For this reason, it can be said it is the most reddish.

It is possible to evaluate differences in brightness and slight differences in color.

## Color Measurement of Eye shadows

Color calculation results (D<sub>65</sub> light source, 2 degree visual field)

No.	Sample name	L*	a*	b*
4	Pink Eye shadow	81.98	15.92	5.33
5	Green Eye shadow	84.22	-12.93	34.01
6	Light blue Eye shadow	79.03	-10.06	-8.57
7	Violet Eye shadow	55.14	6.47	-21.70

The results of color measurement are shown in the L\* a\* b\* colorimetric system, therefore the difference in brightness and color can be easily checked.